

## DuraHand



**DuraHand is a reusable, X-ray compatible, lead free hand immobilizer. The DuraHand material is made of pliable aluminum for easy positioning.**

- **Lead and Latex Free**
- **Allows for the use of X-Ray and Fluoroscopy**
- **Anti-Tilt Tab prevents arm and hand rotation during surgery**
- **Softer and Smoother than the dated lead hand**

Catalog#	Description
<b>DH001</b>	DuraHand: Reusable Hand Immobilizer, Large, 1 Each
<b>DH002</b>	DuraHand: Reusable Hand Immobilizer, Extra Large, 1 Each

P.O. Box 35565 Tulsa, OK 74153-0565

Toll Free 800-948-7994 FAX 918-294-1065

Website: [www.surgicaladvantage.net](http://www.surgicaladvantage.net) Email: [order@surgicaladvantage.net](mailto:order@surgicaladvantage.net)

# Alumi-Hand



**Alumi-Hand is a safe, sterile, ready to use surgical hand immobilizer. The Alumi-Hand is made with a pliable aluminum interior for easy positioning and a foam exterior for patient comfort.**

- **Lead and Latex Free**
- **Allows use with X-Ray and Fluoroscopy**
- **Anti-Tilt Tab prevents arm and hand rotation during surgery**

Catalog#	Description
SA760	Alumi-Hand: Disposable Hand Immobilizer, Large, 8/Box
SA761	Alumi-Hand: Disposable Hand Immobilizer, Extra Large, 8/Box