

...for centuries mariners used Polaris, the North Star, to guide them to their destination. Now, Kirwan Surgical, the world leader in non-stick bipolar forceps, introduces POLARIS. A non-stick bipolar so advanced it resets the standard by which all other non-stick forceps will be measured.

POLARIS™

Your guide to a better surgical experience.



REUSABLE NON-STICK BIPOLAR FORCEPS

POLARIS forceps employ a unique patented mixed matrix metal alloy. This remarkable alloy produces a super conductive bipolar forceps which maintains its properties and non-stick effect at any temperature. The non-stick effect of POLARIS forceps is unlike any bipolar forceps currently on the market and raises the bar on non-stick bipolar forceps performance to a new level.




POLARIS Performance

Kirwan Surgical has been the main innovator of non-stick bipolar technology for forty years. It is the foundation by which Kirwan has maintained its leadership role in the industry. POLARIS forceps represent the next step in Kirwan's non-stick journey. These forceps are extremely conductive, have superior heat-sinking properties while minimizing the residual heat signature in the target tissue. POLARIS forceps will outperform other non-stick forceps and guide you to a better surgical experience.

POLARIS Non-Stick Effectiveness

POLARIS has shown to be more effective in maintaining its non-stick properties in side by side testing against other bipolar forceps. The science is in the material – the unique mixed matrix alloy. In continuous testing, POLARIS forceps maintained a substantially improved non-stick effect and efficiently moved heat away from target tissue when coagulation is achieved.



POLARIS Patterns

The POLARIS reusable forceps are available in four working lengths, 7 inch, 8 inch, 8 ¾ inch, and 10 inch, each with four different tip sizes, 0.5mm, 0.7mm, 1.0mm, and 1.2mm.

Item #	Description
22-1432PL	7" (17.8cm) POLARIS Non-Stick BP Forceps, Ins., 0.5mm
22-1433PL	7" (17.8cm) POLARIS Non-Stick BP Forceps, Ins., 0.7mm
22-1434PL	7" (17.8cm) POLARIS Non-Stick BP Forceps, Ins., 1.0mm
22-1435PL	7" (17.8cm) POLARIS Non-Stick BP Forceps, Ins., 1.2mm
22-1442PL	8" (20.3cm) POLARIS Non-Stick BP Forceps, Ins., 0.5mm
22-1443PL	8" (20.3cm) POLARIS Non-Stick BP Forceps, Ins., 0.7mm
22-1444PL	8" (20.3cm) POLARIS Non-Stick BP Forceps, Ins., 1.0mm
22-1445PL	8" (20.3cm) POLARIS Non-Stick BP Forceps, Ins., 1.2mm
22-1452PL	8 ¾" (22.2cm) POLARIS Non-Stick BP Forceps, Ins., 0.5mm
22-1453PL	8 ¾" (22.2cm) POLARIS Non-Stick BP Forceps, Ins., 0.7mm
22-1454PL	8 ¾" (22.2cm) POLARIS Non-Stick BP Forceps, Ins., 1.0mm
22-1455PL	8 ¾" (22.2cm) POLARIS Non-Stick BP Forceps, Ins., 1.2mm
22-1462PL	10" (26cm) POLARIS Non-Stick BP Forceps, Ins., 0.5mm
22-1463PL	10" (26cm) POLARIS Non-Stick BP Forceps, Ins., 0.7mm
22-1464PL	10" (26cm) POLARIS Non-Stick BP Forceps, Ins., 1.0mm
22-1465PL	10" (26cm) POLARIS Non-Stick BP Forceps, Ins., 1.2mm