

Organize and store your 9" k-wires and pins for easy access and identification.

The gRack holds popular 9" k-wire and pin diameters.

Accommodates 6 each of the following diameters:

0.7mm [.028"]	2.4mm [3/32"]
0.9mm [.035"]	2.8mm [7/64"]
1.1mm [.045"]	3.2mm [1/8"]
1.4mm [.054"]	3.5mm [9/64"]
1.6mm [.062"]	4.0mm [5/32"]
2.0mm [5/64"]	4.5mm [.177"]

- Folds closed for convenient sterilization and storage.
- K-Wires and pins are held in place when rack is closed, helping to prevent shifting of contents.
- Overall size in closed position: 12 1/2"(l) x 5 1/4"(w) x 1 1/2"(h)
- When open, converts to a table top stand for use in the O.R.
- K-Wire and pin diameters are marked on measuring gauge and inside rack for easy identification.
- Easily portable with handle for convenient carry and transport.
- Made from anodized aluminum which is lightweight, yet built to withstand rigorous and repeated use.
- gSource k-wires and pins are available and sold separately. For a complete listing, see the back of this page.

gS 98.5409 holds 9" k-wires and pins
(k-wires and pins sold separately - listed on back)

gRack, K-Wire and Pin

Distributed by:

Surgical Advantage

Phone:
800.948.7994
Fax:
918.294.1065



closed position



open position as a table top stand

Made In Germany.

Unsurpassed quality.

Certified implant steel.

K-Wires

Double Trocar - Smooth						
	diameter	4"	5"	6"	9"	12"
0.7mm [.028"]	78.2000	78.2050	78.1210	78.2105	78.2200	
0.9mm [.035"]	78.2010	78.2060	78.1220	78.2110	78.2210	
1.1mm [.045"]	78.2020	78.2070	78.1230	78.2120	78.2220	
1.4mm [.054"]	78.2040	78.2080	78.1240	78.2140	78.2240	
1.6mm [.062"]	78.2030	78.2090	78.1250	78.2130	78.2230	

Double Trocar - Full Thread

Double Trocar - Full Thread						
	diameter	4"	5"	6"	9"	12"
1.6mm [.062"]	78.4210	78.4220	78.4230	78.4030	78.4035	

Steinmann Pins

Double Trocar - Smooth			Full Thread
	diameter	9" smooth	
2.0mm [5/64"]	78.5500	78.8500	
2.4mm [3/32"]	78.5530	78.8530	
2.8mm [7/64"]	78.5560	78.8560	
3.2mm [1/8"]	78.5590	78.8590	
3.5mm [9/64"]	78.5620	78.8620	
4.0mm [5/32"]	78.5650	78.8650	
4.5mm [1/77"]	78.5680	78.8680	

- Precision ground from certified implant stainless steel.
 - Smooth tapered points are expertly machined for easier penetration.
 - Please inquire about the availability of any size and style not shown on this page.
- Kirschner Wires and Steinmann Pins**
6 per package (non-sterile)

K-Wires

Double Diamond - Smooth						
	diameter	4"	5"	6"	9"	12"
0.7mm [.028"]	78.3000	78.3050	78.1300	78.3100	78.3200	
0.9mm [.035"]	78.3010	78.3060	78.1310	78.3110	78.3210	
1.1mm [.045"]	78.3020	78.3070	78.1320	78.3120	78.3220	
1.4mm [.054"]	78.3030	78.3080	78.1340	78.3140	78.3230	
1.6mm [.062"]	78.3040	78.3090	78.1330	78.3130	78.3240	

Single Diamond - Smooth

Single Diamond - Smooth						
	diameter	4"	5"	6"	9"	12"
0.7mm [.028"]	78.3300	78.3341	78.3350	78.3400	78.3500	
0.9mm [.035"]	78.3310	78.3342	78.3360	78.3410	78.3510	
1.1mm [.045"]	78.3320	78.3344	78.3370	78.3420	78.3520	
1.4mm [.054"]	78.3330	78.3346	78.3380	78.3440	78.3525	
1.6mm [.062"]	78.3340	78.3348	78.3390	78.3430	78.3530	

Steinmann Pins

Double Diamond - Smooth			Full Thread
	diameter	9" smooth	
2.0mm [5/64"]	78.7000	78.8300	
2.4mm [3/32"]	78.7030	78.8330	
2.8mm [7/64"]	78.7060	78.8360	
3.2mm [1/8"]	78.7090	78.8390	
3.5mm [9/64"]	78.7120	78.8420	
4.0mm [5/32"]	78.7150	78.8450	
4.5mm [1/77"]	78.7180	78.8480	

Single Diamond - Smooth

Single Diamond - Smooth			Threaded
	diameter	9" smooth	
2.0mm [5/64"]	78.7780	78.8000	
2.4mm [3/32"]	78.7630	78.8030	
2.8mm [7/64"]	78.7660	78.8060	
3.2mm [1/8"]	78.7690	78.8090	
3.5mm [9/64"]	78.7720	78.8120	
4.0mm [5/32"]	78.7750	78.8150	
4.5mm [1/77"]	78.7782	78.8180	

K-Wires

Single Trocar - Smooth						
	diameter	4"	5"	6"	9"	12"
0.7mm [.028"]	78.2300	78.2700	78.2800	78.2500	78.2600	
0.9mm [.035"]	78.2310	78.2710	78.2810	78.2510	78.2610	
1.1mm [.045"]	78.2320	78.2720	78.2820	78.2520	78.2620	
1.4mm [.054"]	78.2330	78.2740	78.2840	78.2540	78.2640	
1.6mm [.062"]	78.2340	78.2750	78.2850	78.2530	78.2630	

S. Trocar Partial Thread 25mm

S. Trocar Partial Thread 25mm						
	diameter	4"	5"	6"	9"	12"
1.6mm [.062"]	78.9110	78.9112	78.9114	78.9116	78.9118	

S. Trocar Full Thread

S. Trocar Full Thread						
	diameter	4"	5"	6"	9"	12"
1.6mm [.062"]	78.4080	78.4085	78.4088	78.4090	78.4100	78.4110

Steinmann Pins

Single Trocar - Smooth			Threaded
	diameter	9" smooth	
2.0mm [5/64"]	78.6100	78.8700	
2.4mm [3/32"]	78.6130	78.8730	
2.8mm [7/64"]	78.6160	78.8760	
3.2mm [1/8"]	78.6190	78.8780	
3.5mm [9/64"]	78.6220	78.8820	
4.0mm [5/32"]	78.6250	78.8850	
4.5mm [1/77"]	78.6280	78.8880	

gSource part numbers are listed on this page by k-wire or pin style, overall length, diameter. gSource also offers a full line of high quality stainless steel wire and pin management instruments. See Sections 81 through 83 of the gSource catalog on-line at www.gSource.com.