

*stealth
surgical*

**MINIMALLY
INVASIVE**
RETRACTOR SYSTEM



SYSTEM BENEFITS

- Ports & Dilators are 100% titanium with high strength to weight ratio and greater corrosion resistance
- Smaller surgical areas and minimal stress
- Non-magnetic and less reflective in 3D imaging and ideal for intraoperative imagery
- Available in a wide variety of diameters and depths for precise retraction
- Color coded for quick and easy distinction of sizes
- Titanium Flexible Retractor Arm or a Stainless Steel Rigid Retractor Arm offer a variety of table mounts





SS066636
Titanium Flexible Retractor Arm
(Uses SS1671631 Table Post
and SS1671633 Table Clamp)



SS066637
Stainless Steel Rigid Retractor Arm
(Uses SS1671633 Table Clamp)

TITANIUM RETRACTOR SIZES

DIAMETER



LENGTH

3cm

SS1671771

SS1671781

SS1671791

SS1671800

SS1671815

SS1671822

SS1671901

4cm

SS1671772

SS1671782

SS1671792

SS1671801

SS1671816

SS1671840

SS1671902

5cm

SS1671773

SS1671783

SS1671793

SS1671802

SS1671817

SS1671841

SS1671903

6cm

SS1671774

SS1671784

SS1671794

SS1671803

SS1671818

SS1671842

SS1671904

7cm

SS1671775

SS1671785

SS1671795

SS1671804

SS1671819

SS1671843

SS1671905

8cm

SS1671776

SS1671786

SS1671796

SS1671805

SS1671820

SS1671844

SS1671906

9cm

SS1671777

SS1671787

SS1671797

SS1671806

SS1671821

SS1671845

SS1671907

10cm

SS1671920

SS1671930

SS1671940

SS1671950

SS1671960

SS1671970

SS1671980

11cm

SS1671921

SS1671931

SS1671941

SS1671951

SS1671961

SS1671971

SS1671981

12cm

SS1671922

SS1671932

SS1671942

SS1671952

SS1671962

SS1671972

SS1671982

**dilator
required**

SS066620
SS1671531
SS1671532
SS1671537

SS1671538

SS1671539

SS1671540

SS1671541

SS1671542

SS1671909

SS1671797
Standard
Titanium
Retractor



SS1671797B
Beveled
Titanium
Retractor
(add "B" suffix
to standard
part number)



SYSTEM INSTRUMENTATION

All instrument components come standard with Panther ceramic coating.



IVD RONGEURS

Working length 7" (180mm)

Item	Description	
SS1671676E	2mm Micropituitary	
SS1671669E	2mm Straight Pituitary	
SS1671675E	2mm Upbiting Pituitary	
SS1671647E	2mm Downbiting Micropituitary	
SS1671678E	2mm Upbiting Micropituitary	
SS1671670E	2mm Downbiting Pituitary	
SS1671677E	2mm Pituitary With Tooth	
SS1671647E	4mm Straight Pituitary	
SS1671681E*	4mm Upbiting Pituitary	



BAYONET SPURLING KERRISONS

Working length 6" (150mm)

Item	Description	
SS64421	40° Up, 150mm wl, 1mm, non-ejector	
SS64022	40° Up, 150mm wl, 2mm, ejector	
SS64023	40° Up, 150mm wl, 3mm, ejector	
SS64024*	40° Up, 150mm wl, 4mm, ejector	
SS69421	90° Up, 150mm wl, 1mm, non-ejector	
SS69022	90° Up, 150mm wl, 2mm, ejector	
SS69023	90° Up, 150mm wl, 3mm, ejector	
SS69024	90° Up, 150mm wl, 4mm, ejector	



CURETTES

Working length 5 1/2" (140mm)

Item	Description	
SS1671733E	Forward Straight, 1.8mm	
SS1671737E	Forward Angled, 1.8mm	
SS1671740E	Reverse Angled, 1.8mm	
SS1671732E	Forward Straight, 3.6mm	
SS1671736E	Forward Angled, 3.6mm	
SS1671730E*	Reverse Angled, 3.6mm	
SS1671731E	Forward Straight, 5.2mm	
SS1671735E	Forward Angled, 5.2mm	
SS1671739E	Reverse Angled, 5.2mm	

ADDITIONAL

Item	Description	
SS1671680E	Straight Micro Scissors	Working length 7" (180mm)
SS1671679E*	Curved Scissors	Working length 7" (180mm)
SS1671666E	Nerve Root Retractor	Working length 6 1/2" (165mm)
SS1671672E	Suction Retractor	Working length 6" (150mm)
SS167167E	Suction Retractor, Extra Wide	Working length 6" (150mm)
SS1671665E	Wide Nerve Root Retractor	Working length 6 1/4" (160mm)
SS1671646E	Suction Tube #8	Working length 5 1/2" (140mm)
SS1671645E	Suction Tube #10	Working length 5 1/2" (140mm)
SS1671635E	Suction tube #12	Working length 5 1/2" (140mm)

SYSTEM INSTRUMENTATION **cont.**

All instrument components come standard with Panther ceramic coating.



90° BALL PROBES

Working length 5" (125mm)

Item	Description	
SS1671752E	Short, Right	
SS1671753E	Short, Straight	
SS1671754E	Short, Left	
SS1671755E	Long, Right	
SS1671756E*	Long, Straight	
SS1671757E	Long, Left	

DISSECTORS

Working length 5" (125mm)


Item	Description	
SS1671758E	Penfield #4 Push	
SS1671759E	Penfield #4 Pull	
SS1671750E	Penfield #2 Push	
SS1671761E	Penfield #2 Pull	
SS1671762E	Nerve Hook, Right	
SS1671763E	Nerve Hook, Straight	
SS1671764E	Nerve Hook, Left	
SS1671765E*	90° Dissector, Right	
SS1671766E	90° Dissector, Straight	
SS1671767E	90° Dissector, Left	
SS1671751E	Woodson Probe 45°	

TRAYS


The Minimally Invasive Port System can be customized to meet your needs. Select the retractor sizes, instruments and trays most relevant to your procedures or work with your local Stealth representative to create the trays best suited for your practice.

Pre-Built Trays


Item	Description
SST2930	MIS Retractor Case I, 3-level, incl 73-2931, 73-2932, 73-2933
SST2935	MIS Retractor Case II, 4-level, incl 73-2936, 73-2937, 73-2938 (x2)
SST2939	MIS Retractor Case For 73-2956 Rigid Arm And 73-2955 Table Clamp
SST2931*	MIS Retractor Tray, Dilators and Ports
SST2932	MIS Retractor Tray, Curettes and Larger Dilators
SST2933*	MIS Retractor Tray, IVD Rongeurs
SST2936	MIS Retractor Tray, Bayonet Dissectors
SST2937	MIS Retractor Tray, Suctions and Retractors
SST2938	MIS Retractor Tray, Holds 3 Bayonet Kerrisons



SST2931
MIS Retractor Tray
Dilators and Ports



SST2933
MIS Retractor Tray
IVD Rongeurs



***stealth
surgical***

104 Sommerfield Drive
Gordonsville, VA 22942 USA

Toll Free: 866-473-2004
Phone: 540-832-5580
Fax: 540-832-5585

www.stealthsurgical.com

©2020 Stealth Surgical, LLC REV B L83
Stealth and Stealth Surgical LLC are Registered Trademarks