




GelPro

Item Code	Size	Color
SA1824GB	18 x 24	Columbia Blue
SA2032GB	20 x 32	Columbia Blue
SA2048GB	20 x 48	Columbia Blue



Item Code	Size	Color
SA1824GR	18 x 24	DND Red
SA2032GR	20 x 32	DND Red
SA2048GR	20 x 48	DND Red



- ¾" Dual-Density™ comfort core featuring proprietary gel + energy-return foam
- Anti-microbial and anti-bacterial properties for added protection
- Water-resistant and will not absorb fluids
- Proper friction level for damp or dry surgical booties
- Easy-to-clean and compatible with hospital grade quaternary cleaners
- NFPA 260 and CAL 117 compliant
- Certified by the National Floor Safety Institute for the highest level of traction
- Available in DND Red for use in no-interruption zones such as medicine dispensing stations
- Made in USA

P.O. Box 35565 Tulsa, OK 74153-0565

Toll Free 800-948-7994 FAX 918-294-1065

Website: www.surgicaladvantage.net Email: order@surgicaladvantage.net

Eco-Pro

Item Code	Size	Color
SA2032EK	20 x 32	Black
SA2048EK	20 x 48	Black
SA2436EK	24 x 36	Black
SA2072EK	20 x 72	Black
SA3660EK	36 x 60	Black

Item Code	Size	Color
SA2032ET	20 x 32	Taupe
SA2048ET	20 x 48	Taupe
SA2436ET	24 x 36	Taupe
SA2072ET	20 x 72	Taupe
SA3660ET	36 x 60	Taupe



Beveled Edge
beveled at 18° to reduce tripping potential


Matte Textured
top surface provides optimum friction

Durable
commercial grade heavy-duty polyurethane

Anti-microbial
properties for added protection

Eco-Friendly
constructed with Bio-Foam® a renewable, plant-based resource

3/4" thick

NFSI  **High-Traction Bottom**
certified by the National Floor Safety Institute

Integral Skin
stain & puncture-resistant surface prevents absorption of dirt & water

Energy-Return Comfort
proprietary Cellulon® core technology with energy-return properties

P.O. Box 35565 Tulsa, OK 74153-0565

Toll Free 800-948-7994 FAX 918-294-1065

Website: www.surgicaladvantage.net Email: order@surgicaladvantage.net