

Integra® Jarit®

Sterilization Containers

Limit uncertainty by protecting
your instrument investment



INTEGRA®
LIMIT UNCERTAINTY

Summary of Integra Sterilization Container Product Features

Gaskets within retention plates and container perimeter ensure a secure seal of the container, lid and filter.

Snap on filter plate with a simple locking mechanism for quick, secure and easy change of disposable filters.

Round disposable filters with centering hole to assist in proper alignment.

Lid and Bottom manufactured from anodized aluminum alloy. Easy to maintain, corrosion resistant, lightweight yet durable, and conducts heat for material drying period.

4 ID Tag Windows (2 front & 2 back)

ID tags available in 6 colors and can be engraved with up to 18 characters per line.

Low profile latch easily and securely closes lid over container bottom.

Tamper Evident Lock – identifies sterility breach.

Fold-down handles extend and stop at 90° for a secure grasp to aid lifting and avoid pinching.

Corner design allows for same size container stacking in storage.

2 Indicator Card Windows (1 front & 1 back). Indicator card has visible bar for clear indication that steam or gas sterilization has been achieved.

*Baskets are constructed with solid wire mesh to allow for maximum ultrasonic cleaning; handles can be positioned to assist with a secure grasp for removal under sterile conditions.

One-Half Size (11¼" x 11") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4"	730-400	730-400B	730-410	730-410B
5⅓"	730-401	730-401B	730-411	730-411B
6"	730-402	730-402B	730-412	730-412B
8"	730-403	730-403B	730-413	730-413B
10"	730-404	730-404B	730-414	730-414B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



One-Half Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-120
Blue	730-121
Red	730-122
Gold	730-123
Green	730-124
Black	730-125



One-Half Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-400	750-400
3"	740-401	750-401
4"	740-402	750-402
Lid	-	750-420



Wire Basket



Perforated Basket

Three-Quarter Size (18¼" x 11") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4"	730-500	730-500B	730-510	730-510B
5⅓"	730-501	730-501B	730-511	730-511B
6"	730-502	730-502B	730-512	730-512B



Three-Quarter Size Colored Lids
(bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-220
Blue	730-221
Red	730-222
Gold	730-223
Green	730-224
Black	730-225




Three-Quarter Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-500	750-500
3"	740-501	750-501
4"	740-502	750-502
Lid	-	750-520



Wire Basket



Perforated Basket

Full-Size (23" x 11") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4"	730-600	730-600B	730-610	730-610B
5 1/3"	730-601	730-601B	730-611	730-611B
6"	730-602	730-602B	730-612	730-612B
8"	730-603	730-603B	730-613	730-613B
10"	730-604	730-604B	730-614	730-614B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



Full-Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-320
Blue	730-321
Red	730-322
Gold	730-323
Green	730-324
Black	730-325



Full-Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-600	750-600
3"	740-601	750-601
4"	740-602	750-602
Lid	-	750-620



Wire Basket



Perforated Basket

Wide-Body (23⅓" x 15¾") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
120mm	730-701	730-701B	730-711	730-711B
180mm	730-702	730-702B	730-712	730-712B



Wide-Body Lid (bottoms sold separately)


Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-220



Wide-Body Wire Baskets


Height	Wire Basket
2"	740-700
4"	740-701



Instrument Stringers and Holders

Catalog Number	Description	Size
117-240	Weinstein Instrument Holder	6"
117-242	Weinstein Instrument Holder	8"
117-244	Weinstein Instrument Holder	10"
117-246	Weinstein Instrument Holder	12"
117-247	Weinstein Instrument Holder	14"
117-248	Weinstein Instrument Holder	16"
117-249	Weinstein Instrument Holder	18"
117-250	Jarit Instrument Holder	6"
117-252	Jarit Instrument Holder	8"
117-254	Jarit Instrument Holder	10"
117-256	Jarit Instrument Holder	12"
117-258	Jarit Instrument Holder	14"
117-260	Jarit Instrument Holder	16"
117-262	Jarit Instrument Holder	18"
117-270	Instrument Stringer	8"
117-272	Instrument Stringer	10"
117-274	Instrument Stringer	12"
117-276	Instrument Stringer	14"
117-278	Instrument Stringer	16"
117-280	Instrument Stringer	18"

Note: For more instruments stringers please refer to the Jarit® catalog.



Weinstein Instrument Holder

Jarit Instrument Holder

Instrument Stringer

Identification Tags

Each tag allows for 2 lines of custom text, 18 characters per line.

Catalog Number	Color
775-100	Silver
775-101	Blue
775-102	Red
775-103	Gold
775-104	Green
775-105	Black



Snap-on Basket Tags (with key ring)

Each small tag allows for 3 lines of custom text, 32 characters per line.
Each large tag allows for 5 lines of custom text, 32 characters per line.

Color	Small	Size	Large	Size
Silver	775-120	2 3/8" x 1/2"	775-130	2 3/8" x 1 1/4"
Gold	775-121	2 3/8" x 1/2"	775-131	2 3/8" x 1 1/4"
Red	775-122	2 3/8" x 1/2"	775-132	2 3/8" x 1 1/4"
Blue	775-123	2 3/8" x 1/2"	775-133	2 3/8" x 1 1/4"
Green	775-124	2 3/8" x 1/2"	775-134	2 3/8" x 1 1/4"
Black	775-125	2 3/8" x 1/2"	775-135	2 3/8" x 1 1/4"




Tamper Evident Locks and Indicator Cards

Catalog Number	Description
775-400	Tamper Evident Locks, single use, blue, box of 1000
775-300	Indicator Cards, box of 250




Tamper Evident Lock




Indicator Card

Retention Plate and Filters

Catalog Number	Description
775-212	Retention Plate
775-200	Disposable paper filters, single use, 7 1/2", round, box of 1000, GRAVITY
775-250	Disposable paper filters, single use, 7 1/2", round, box of 1000, PRE-VACUUM



Retention Plate




Disposable Paper Filters

Instrument Rack

- 9 slots
- 5-10mm diameter silicone bars
- Compatible with the 8" Full-Size container
- The rack fits inside the 2" or 3" wire or perforated basket, allowing ample storage underneath for other instruments and auxiliary items

Catalog Number	Description
790-010	Instrument Rack



Silicone Holders

Catalog Number	Description	Size
775-630	Silicone Holder #1	8mm x 26mm x 40mm
775-631	Silicone Holder #2	8mm x 26mm x 40mm
775-632	Silicone Holder #3	8mm x 26mm x 40mm
775-633	Silicone Holder #4	10mm
775-634	Silicone Holder #5	5mm
775-650	Silicone Holder Clamps, 4 slot, plastic	



#1



#2



#3




#4



#5



Silicone Mats

Catalog Number	Description
776-400	One-Half Size Silicone Mat
776-500	Three-Quarter Size Silicone Mat
776-600	Full-Size Silicone Mat



Partition Sheets and Holding Pins

Catalog Number	Description	Size
775-610	Partition Sheet	2" x 7/8" (50mm x 20mm)
775-611	Partition Sheet	5" x 7/8" (130mm x 20mm)
775-612	Partition Sheet	9 7/8" x 7/8" (225mm x 20mm)
775-613	Partition Sheet	18 1/8" x 7/8" (460mm x 20mm)
775-600	Holding Pin for Wire Baskets	
775-601	Holding Pin for Partition Sheets	




Holding Pin for Wire Baskets



Holding Pin for Partition Sheets

Holding Clamps

Catalog Number	Size Range
775-620	1/8" to 1/2" (7mm-12mm)
775-621	1/2" to 5/8" (12mm-16mm)
775-622	5/8" to 1" (16mm-26mm)
775-623	1" to 1 1/2" (26mm-36mm)
775-624	1 1/2" to 1 3/4" (36mm-45mm)



Tolerances in product specifications, due to the hand craftsmanship of the instruments, are subject to minor variations in dimensions and finishes.

For more information or to place an order, please contact:

Integra ■ 589 Davies Drive, York, PA 17402

877-468-5572 USA ■ 914-789-7094 outside USA

609-750-4259 fax ■ 914-592-8056 outside USA fax

integralife.com

