

- Instrument stringer holds finger ring handle instruments securely ensuring instruments stay open during sterilization.
- Locking ball closure keeps cross bar in place.
- Stay-closed design prevents accidental opening.
- Expertly hand finished to eliminate any sharp edges.
- Made from German stainless steel.
- To open or close: squeeze side bars together, then fold back cross bar.




	2 1/2" wide		3 1/2" wide
<b>gS 98.2104</b>	4"		
<b>gS 98.2106</b>	6"	<b>gS 98.2206</b>	6"
<b>gS 98.2108</b>	8"	<b>gS 98.2208</b>	8"
<b>gS 98.2110</b>	10"		
<b>gS 98.2112</b>	12"		
<b>gS 98.2114</b>	14"		

**Instrument Stringer**  
with lock

<b>gS 98.2308</b>	8"
<b>gS 98.2310</b>	10"
<b>gS 98.2312</b>	12"
<b>gS 98.2314</b>	14"
<b>gS 98.2316</b>	16"
<b>gS 98.2318</b>	18"

**Instrument Stringer**  
2 3/4" wide  
U-shaped, open end



<b>gS 98.2406</b>	6"
<b>gS 98.2408</b>	8"
<b>gS 98.2410</b>	10"
<b>gS 98.2412</b>	12"

**Weinstein Instrument Holder**  
2 3/4" wide  
U-shaped, closed end

