

SKIN TRANSPLANTATION

We care



Technical details Dermatome D80

D80 Dermatome Features

- Cordless, lightweight dermatome for optimal free movement
- Easy thickness adjustment from 0.0 to 1.2mm with 0.1mm increments
- Standard cutting width of 80mm, supplied with reducing clamps for a width of 35, 50 and 65 mm
- Quick blade replacement possible without tooling
- Anodised aluminium and stainless steel

Technical data

Weight	1330 g / 47 oz
Dimensions (lxw)	295 x 104 mm / 11.6" x 4.09"
Diameter grip	45 mm / 1.77"
Cutting width	80 mm / 3.15"
Graft thickness	0-1.2 mm / 0-0.048"
Stroke blade	3 mm / 0.118"
Dimensions sterilization case (lxwxh)	480 x 250 x 66 mm / 18.9" x 9.8" x 2.6"
Motor capacity	15 W
Battery (Li-Ion)	7.4 V / 2400 mAh
Charge supply	100-240 V
Charge current	1200 mA
Charge time	~ 3 hrs

Set contents

4.D80STS	D80 dermatome, complete set, 2400mAh Li-Ion battery, this includes:
4.D80	D80 dermatome without battery- and motor cartridge
4.ANS-HC	Dermatome battery charger with adaptors
4.SF02	Sterile funnel for large diameter dermatome shaft
4.SCL	Sterile clamp
4.D80AC2	D80 sterilization case
4.D80CL35-20	D80 width reducing clamp 35 mm
4.D80CL50-21	D80 width reducing clamp 50 mm
4.D80CL65-22	D80 width reducing clamp 65 mm
4.SU01	Dermatome charger support unit
4.BC7.4V2400	Dermatome battery cartridge Li-Ion, 7.4V, 2400 mAh
4.MCX15	Dermatome motor cartridge, long (new) version
4.076	Torx T20 screwdriver
5.D80BL10	D80 dermatome blades, box 10 pcs.



4.D80STS



4.D80



4.ANS-HC



4.SF02




4.SCL



4.D80AC2



4.D80CL35 / 4.D80CL50 4.D80CL65



4.SU01



4.BC7.4V2400



4.MCX15



5.D80BL10



4.076