

POLARIS[®]

Bipolar Electrosurgical Generator

The **POLARIS** bipolar generator by Kirwan Surgical Products represents the newest technology available in the marketplace today. This state-of-the-art electrosurgical generator with modular accessories boasts a fully interactive touch screen with drop-down menus for programmable user profiles, providing the ultimate in user customization and ease of use.



29-1000 POLARIS Bipolar Generator

The **POLARIS** bipolar generator is field upgradable, ensuring your equipment stays up-to-date with the latest advancements in surgical technology. Its flexible modular design allows for complete user control of all functions, including bipolar coagulation, bipolar cut, irrigation, and illumination.

- Modular design capable of supporting up to four independent modules
- LCD capacitive touch screen stores up to 31 profiles and allows users to quickly adjust and fine-tune generator settings
- Continuity bipolar accessory indicator ensures proper connection
- Separate output settings for both cut and coagulation waveforms
 - Coagulation sinusoidal waveform with specific duty cycle for more efficient tissue coagulation
 - Cut sinusoidal waveform for a more effective outcome



Order Information:

Product #	Description
29-1000	Polaris Generator
29-1600	Polaris Irrigation Module
29-1900	Polaris Light Module
10-4603	Polaris Irrigating Tubing Set

Also available:

Disposable Polaris Forceps

- 3 lengths
- 3 tip sizes



Reusable Polaris Forceps

- 4 lengths
- 4 tip sizes

